

First Class Laser
Made in USA

Cellulaze®



CYNOSURE®

Amico Yasna Pars (Pr. J. S. Co)

Cellulaze®

Nd:YAG Laser

CYNOSURE®
A Hologic Company

The powerful Nd:YAG laser provides selective and lasting treatment of the structural causes of cellulite in deep layers of tissue

INDICATIONS

- **Cellulite**
- **Body Contouring**

PRODUCT ADVANTAGES

- **Precise and long lasting results**
- **Minimally invasive**
- **Fights the three causes of cellulite under the skin:**
Smooths pads of fat, thermally separates fibrous septae and firms the tissue by stimulating the production of collagen and elastin crossing the hypodermis
- **Unique SideLaze800™ fibre with laser beam in 2 directions**

Cellulite, skin tightening



Picture source: Dr. G. Sasaki



Picture source: Dr. G. Sasaki

TECHNICAL DATA

System	Nd:YAG
Wavelength	1440 nm
Pulse duration (max.)	1000 µm
Output	10 W
Max. repetition rate	25 Hz
Laser beam	red
Electric requirement	230 VAC ± 10%, 16 A
Size (H x W x D)	95 cm x 50 cm x 65 cm
Weight	100 kg