

First Class Laser
Made in USA

PicoSure®



Focus™



CYNOSURE®

Amico Yasna Pars (Pr.J.S.Co)

PicoSure® is the only 755 nm/532 nm/1064 nm laser for aesthetic skin treatments and tattoo removal of all colors. PicoSure® Focus™ uniquely treats pigmented lesions, acne scars and wrinkles with impressive results without the typical downtime or discomfort. PicoSure® is the world's first picosecond aesthetic laser platform widely recognized for its technology leadership, clinical versatility and proven performance. PicoSure® has enabled clinical applications not available with any other laser

INDICATIONS

- **Tattoo removal*** of any color
- **Pigmented lesions***
- **Acne scars***
- **Wrinkles***
- **Skin rejuvenation**

(*FDA approved)


PRODUCT ADVANTAGES

- **Fewer treatments, faster clearance**
- **25 - 50% fewer sessions required vs. Q-Switched lasers**
- **Photomechanical efficacy, fewer side effects**
- **Proven performance across several indications**

TECHNICAL DATA


System	Alexandrite/picosecond pulsed
Wavelength	755 nm (standard) 532 nm (optional) 1064 nm (optional)
Pulse duration	550 - 750 ps
Pulse power	165 - 200 mJ
Handpieces	Zoom with continuous adjustment from 2.0 to 6.0 mm Additional fixed handpieces 6, 8 and 10 mm
Max. repetition rate	10 Hz
Pulse duration	550 - 750 ps
Electrical power	380 V~, 4,5 kVA, 50/60 Hz, single-phase
Electrical requirement	30 A
Size (H x W x D)	56 cm x 107 cm x 158 cm
Weight	171 kg

Tattoo removal*



Picture source: Dr. J. Dover und Dr. K. Arndt

Tattoo removal*



Picture source: Zieker Ophthalmology

Acne scars*



Picture source: Dr. R. Saluja

Tattoo removal*




Picture source: Dr. K. Arndt

Pigmented lesions*



Picture source: Zieker Ophthalmology

Wrinkles*



Picture source: Dr. R. Weiss

