

First Class Laser
Made in USA

SculpSure®



WarmSculpting™
with SculpSure®



CYNOSURE®

Amico Yasna Pars (Pr. J. S. Co)

The SculpSure® device provides WarmSculpting™ treatments for **non-invasive body contouring** that permanently reduces stubborn fat without surgery or downtime. WarmSculpting™ for the body (flanks, abdomen, thighs, and back) is intended for clients with a BMI up to 30. Patients contouring the submental area can have a BMI of up to 43

INDICATIONS

- **Fat reduction: abdomen, flanks, back, inner and outer thighs, submental**

PRODUCT ADVANTAGES

- **Fast, 25-minutes treatment per area**
- **Versatile applicators to fit a variety of body shapes and sizes**
- **No surgery, no downtime**
- **High ROI to increase your customers revenues faster**



TECHNICAL DATA

System	Diode Laser
Wavelength	1060 nm
Pulse width	CW
Power output	40 W max per diode (160 W total)
Aiming beam	none
Electric requirement	230 VAC (20 amp circuit required)
Size (HxBxT)	144.8 cm x 61 cm x 91.5 cm" with arm, no umbilical loop 172.8 cm x 61 cm x 91.5 cm with arm and umbilical loop
Weight	118 kg with applicators 112 kg without applicators

