

First Class Laser
Made in USA

Smartlipo[®] Triplex[®]



CYNOSURE[®]

Amico Yasna Pars (Pr. J. S. Co)

Smartlipo® Triplex®

Nd:YAG Laser

The Smartlipo® Triplex® is the “gold standard” of laser lipolysis. It is the only system featuring 3 wavelengths - 1064 nm, 1320 nm and 1440 nm. These wavelengths can be used individually or combined for so-called MultiPlex™ technology. Optimal results are achieved together with the unique SmartSense™ (movement control) and ThermaGuide™ (temperature control)

INDICATIONS

- **Controlled liposuction assisted by laser**

PRODUCT ADVANTAGES

By selecting wavelengths according to the condition and patient, Smartlipo® Triplex® offers:

- **Tissue tightening using soft tissue coagulation**
- **Reliable, controlled energy output**
- **A shorter treatment duration and patient downtime**
- **Less bleeding and haematomas than classic liposuction**
- **Defined body sculpting**
- **Minimally invasive procedure**
- **Constant temperature monitoring**

Body contouring, liposuction and skin tightening



Picture source: Dr. C. Petti



Picture source: Dr. C. Petti

TECHNICAL DATA

System	Nd:YAG laser		
Wavelength	1064 nm	1320 nm	1440 nm
Pulse duration (max.)	500 μs	1000 μs	1000 μs
Output	40 W	24 W	15 W
Max. repetition rate	40 Hz	40 Hz	25 Hz
MultiPlex options	1064 nm and 1320 nm 500 μs and 1000 μs	1064 nm and 1440 nm 500 μs	
		and 1000 μs	
		40 Hz and 25Hz	
		40Hz	
		red	
Laser beam	red		
Electric requirement	230 VAC, 16 A		
Size (H x W x D)	95 cm x 50 cm x 65 cm		
Weight	100 kg		