

If You Wake Up Without Goals, Go Back To Sleep...



Inside this issue:

TECNIS1: Optical Synergy 2

ALCON VISION CARE: 2
Dailies brand and all day comfort technology

HUVITZ: 3
Digital Refractor (Phoropter)

OASIS: New Iris Expander 3

ABBOTT MEDICAL OPTICS 4

24th Annual Congress of the Iranian Society of Ophthalmology 4

This is the motto that we, as **Amico Yasna Pars** staff refer to, everyday.

We all have goals in our life which varies according to our conditions, your goal may be your patients' satisfaction, and our goal is helping you reach there easier, we are walking on same line, but still you may want to ask us that what's the meaning of the motto.

The answer is very simple; Seeing you, Being with you, and using your comments and suggestions as our guide light, is our goal.

Sharing what we have with you, as our partners, and acquiring your knowledge to

improve our technical skills in introducing our products more effective, is the energy we need to continue, in the way of creating new lines for **Amico Iran**, and also proffering better solutions for your patients.

We are always seeking your warm and kind suggestions in order to atone our situation with what is your demand, and we have proven our flexibility in amending the situation to change an impasse to a great highway of communication and relation.

As always, **Amico Team** embrace this moment to congratulate you coming the new year, and also wish you its best wishes.

We are grateful of having you

TECNIS® 1

Abbott
Medical Optics

TECNIS, This is Optical Synergy

Obviously you are worried about your patient's post-op outcomes. One of the main complaints of patient's is having CAT EYE phenomenon in their eyes.

What is the reason?

The reason is presence of micro vacuoles inside the lens material which will reflect light and cause CAT EYE phenomenon.

Also Glistening can cause loss in contrast sensitivity & Decreased visual acuity.


TECNIS IOL's are glistening free, because:

1. Manufactured by lathe cutting which results in fewer gaps vs. molding.
2. **TECNIS®** IOLs are made with a proprietary


cryo-lathing method that limits microvoid formation and high temperature fluctuations for reduced glistening formation

3. Temperature fluctuations of other IOLs have been shown to cause glistening formation.

Acrysof



TECNIS 1 Piece



Thinking about contact lenses?

Alcon
Vision Care

Your patient may be surprised at how comfortable, convenient and affordable contact lenses can be. Today there are a wide range of options for contact lenses so you can give your patient exactly what they need

ALL DAY COMFORT

- Refreshing all day comfort with every blink

All **DAILIES®** brand contact lenses provide moisture with every blink, which is up to 14,000 times per day on average.

ADVANCED TECHNOLOGY

Only **DAILIES®** brand contact lenses moisturize with every blink, right up to the end of day. **DAILIES®** brand contact

lenses fit comfortably into user life, too. Just put on a fresh, new pair each morning and dispose them at night! **DAILIES®** brand contact lenses are all hassle free with no lens care solutions required. These contact lenses are designed to be convenient for users.

Plus all, **DAILIES®** brand contact lenses are made with an incredibly thin edge for extra comfort when user put them on and during wear.

BENEFITS

Freedom with a new lens feeling every day

- Contacts match everything user wear
- Lenses don't slip down user nose
- They don't fog up or get splashed

with rain

- No annoying obstructions or reflections in view
- Great choice for sports
- Can be worn under protective eyewear



Comes with ALCON Packaging...

Huvitz



Digital Refractor (Phoropter)

The **HUVITZ Digital Refractor HDR-7000** system is a computerized auto refractor used for subjective refraction regarding the patient's accommodation power and facilities. The subjective refraction can be practiced on the basis of objective refraction that is done by using an **Auto Ref/ Keratometer (RK)** connected to Digital Refractor.

The examiner can do optometric practice in pretty short period of time with ease with the aid of fully integrated and automated Digital Refraction System, where every instrument for the refraction is interfaced and controlled by one touch of button on the Digital Refractor.

HDR-7000 system includes a refractor body, an operation panel and a junction box. The refractor body has more than a hundred lenses and serves as an automatic lens-loading unit. It is connected to the operation panel via the junction box and is controlled by the operation panel. The operation panel is the main controller of **HDR-7000** system. It controls the refractor body, chart-presenting device and also receives data transmitted from **Auto Ref/Keratometer(RK)**, **Auto Lensmeter(LM)** or a PC interfaced with Digital Refraction System and also built-in thermal printer embedded.

The junction box serves as a power supply to the Digital Refractor and the Operation Panel and it provides communication ports for interfaced instruments such as **Chart Projector** or **Auto Ref/Keratometer**.

Oasis: Looking Forward Surgeons: Look For Changes

OASIS Medical, Inc. is a specialty manufacturer of EYE care products and surgical supplies, instruments and disposables. The OASIS Company founded in 1987 by Norman Delgado and headquartered near Los Angeles.


New Iris Expander Available from Oasis Medical:

Oasis Medical Inc. recently introduced its new 6.25 mm Oasis Iris Expander, which is used for expanding insufficiently dilated pupils during ophthalmic surgery and sustaining visibility throughout the procedure. Once inserted, the molded

polypropylene ring has four pockets that cradle the delicate ocular tissue without pinching or clamping for an unimpeded view during the surgical procedure and easy removal without potential tearing. The Iris Expander is packaged sterile with an inserter and comes with five single-use units per box.

The pupil expander devices currently on the market differ with regard to ease of handling and clinical results. The Iris Expander is a useful tool for

phacoemulsification surgery. It adequately dilates the pupil and prevents iris sphincter damage. This easy-to-insert-and-remove device expands the pupil, protects the iris sphincter during surgery, and allows the pupil to return to its normal shape, size and function after the operation.



OASIS®

Abbott Medical Optics



There are usually some problems with the way you prefer to implant the IOL's. Except the surgeons' preference, **Abbott Medical Optics** has come up with a new solution for **TECNIS** IOL's.

UNFOLDER Platinum Implantation system;

Designed for Micro-Implantation procedures (2.2 - 2.4 mm Incisions) and compatible across the full range of **TECNIS® 1-PIECE** IOL DIOPTERS (5-34D)

As you see, you can use this method for all diopters of all **TECNIS 1 Family of IOL's***

There are different types of Injectors produced by AMO for different tastes:

1. DK 7796, Unfolder Platinum 1 Series Implantation Systems to cover Micro Incision surgeries. With its screw style, you will have a complete control over your implantation.
2. DK 7786, ONE SERIES Ultra Implantation systems with plunger style (one handed) for 2.8 mm and 3.2 mm incisions.
3. EMERALD Screw Style series implantation systems for 3 Piece IOL's

* TECNIS 1 Family of IOL's:

1. TECNIS 1 Monofocal IOL's
2. TECNIS Multifocal IOL's
3. TECNIS Toric IOL's
4. TECNIS Multifocal Toric IOL's

24th Annual Congress of the Iranian Society

The 24th annual congress of the Iranian society of ophthalmology was held in Tehran, Iran, December 1 – 4 at Razi convention center.

Amico Yasna Pars was honored to be one of the main sponsors of the congress and provided a wide range of leading brands products and equipment in the field of ophthalmology and optometry.

We're also pleased to announce the presence of **Abbott Medical Optics** cataract products for the first time in this congress, which attracted ophthalmology society of Iran attention.

During the visit of Dr. Hassan Hashemi, administrator of health ministry of Iran from the congress booths, it was a great chance for different partners in exhibition to introduce their products, and also we have faced with his warm welcomes and wise comments.

Respectfully, **Amico Yasna Pars** invites you to visit its website for reviewing congress photos, ordering marketing materials and to share ideas about congress.
www.amicoyasnarpars.com



Amico Yasna Pars (Pr.J.S.Co)

6th floor, No.1643, Shariati Ave., Tehran, Iran
Postal Code: 1939613881 Tel: (+98)21-22645870-71 Fax: (+98)21-22645872
Email: amicoyasnarpars@amicogroup.com
Website: www.amicoyasnarpars.com